
# KERN **572-55**







| Measuring system                                |  |
|---|--|
| Adjustment options:                             | External calibration   |
| Linearity:                                      | 250 mg   |
| Readability [d]:                                | 50 mg  |
| Recommended adjusting weight:                   | 20 kg (F1)   |
| Repeatability:                                  | 100 mg   |
| Resolution:                                     | 400000   |
| Stabilisation time under laboratory conditions: | 3 s  |
| Warm up time:                                   | 4 h  |
| Weighing capacity [Max]:                        | 20 kg  |
| Weighing system:                                | Strain gauge   |
| Weighing units:                                 | t ozt kg to tl (HK) g mo tl<br>(Cn) lb tl (Singap. Malays)<br>dwt ct oz tl (Tw) gn |
| Display   |  |

| Біэріаў   |          |  |
|---|----------|--|
| Display digit height (large):                       | 1,800 cm |  |
| Counting  |          |  |
| Counting resolution:                                | 400000   |  |
| Mininum piece weight at piece counting (Laboratory: | 50 mg    |  |
| Mininum piece weight at piece counting (Normal):    | 500 mg   |  |
|   |          |  |

| I |
|---|
|   |
|   |
|   |
| l |
|   |
|   |
|   |
|   |
|   |
|   |
|   |
|   |
|   |
|   |
|   |
|   |
|   |
|   |
|   |
|   |
|   |

Revolving screw feet: Weighing surface (WxD):

160 x 200 mm

## KERN **572-55**



All-rounder e.g. as precision balance in the laboratory or in harsh industrial applications

### Piktogramme

### STANDARD



#### **OPTION**

